

Bubble CPAP System

The respiratory humidifier ensures accurate flow with leakage prevention, supports invasive and non-invasive modes with auto temperature control, works with all ventilators, offers pediatric/adult connectors, supports reusable/disposable circuits, includes a built-in alarm, and comes with a safe mounting clamp and slot plate.

- Flow meter from 0.1 to 1LPM and from 1 to 10LPM
- Support dual output with double flow meters
- Made with high quality aluminum alloy, anodic oxidation surface treatment
- Adjust oxygen concentration level from 21% to 100%
- Powered by pneumatic control, no electricity required
- Mechanical alarm system by detection of pressure differential



Air Oxygen Blander



Respiratory Humidifier



Specifications

FiO2	: 21%-100%
Flow Control	: 0.1-1 & 1-10LPM
Pressure	: Air/Oxygen@0.3-0.4MPa
Alarm	: When gas supply pressure difference 0.1MPa
Number of Output	: Dual
Accuracy of FiO2	: $\pm 3\%$
Accessories	: Heated Plate
Power Source	: AC 220V, 50-60Hz :150W
Thermal Cut off	: 95 ± 6
Display	: Three digit, 7 segment Led
Display Size	: 0.56 inch
Temperature Monitoring	: 0-75

Accessories



Bubble Generator



Head Wire Circuit



CPAP Tubing



Infant Bannot



Nasal Prong

GuangDong Pigeon Medical Apparatus Co., Ltd

Add.: No.118 Zhuyuan, Chiwei Town, Punning City, Guangdong, China

Tel: +86(653)2298001

Email: service@gdgezi.com

Website: www.gdgezi.com