

MODEL NO. HB-M202



Manual Two Function Hospital Bed

Configuration Details:

Detachable ABS Head & Foot Board.

Overall Length : 2080mm.

Overall Width : 960mm.

Sleep Surface length : 1960mm.

Sleep Surface Width : 900mm.

Sleep Surface height : 500mm.

Collapsible SS. Guide Rails : 6 Bars Metal joint

Weight of bed : 50Kg.

Load capacity of bed : 240Kg.

Noise-less 5" Castors with Brake.

Manufacturer of Hiking Medical Co. CHINA.