

SV300 Ventilator

Compact yet powerful







www.mindray.com





Design based on simple facts

Comprehensive

SV300 comes with extensive ventilation modes and is ideal for acute ICU use.

Adaptive

With useful functions, it is adaptable to different clinical environments.

Simple

Intuitive user interface optimizes clinical workflow.



We listen to **your opinions** and try best to make your **daily work easier**



Don't be fooled by the compact size of the SV300 ventilator. It comes with extensive ventilation modes and is equipped with functions that are usually only found on intensive care ventilators. We have also integrated unique features that you won't find on other high spec models.



ICU package options

- APRV
- PRVC
- ATRC
- Static PV Loop
- Recruitment tool SI



Extensive monitoring

- 4 waveforms
- 2 loops
- 72 hrs patient trend data
- 5000 events log



Volumetric CO₂ measurement

- Airway dead space
- CO₂ elimination



Weaning indicators

- P0.1
- NIF
- RSBI



O₂ Therapy

- O₂ enrichment
- Adjustable flow



IntelliCycle™ Pro

- Automatically adjusts trigger level, Tslope and ETS
- Improves patient-ventilator synchrony



Versatile functions all on one device



Invasive therapy: Bedside intubated

- Extensive ventilation modes
- ATRC
- Weaning indicators



Non-invasive and O₂ therapy

- Exceptional leak compensation capability
- Quiet turbine
- Integrated SpO₂ monitoring



On the run: Intra-hospital transport

- Easily detachable trolley
- Approximate 10kg
- Up to 6 hours internal battery support



CPRV™ ventilation

- Integrates unique e-ITD™
- Improves venous return and helps improve perfusion
- Automatically associates CO₂ monitoring



Neonatal ventilation

- Comprehensive ventilation modes for neonates
- Accurate proximal flow sensor
- High precision flow and pressure control



- Adjustable screen allows best viewing angle



- PulmoSight[™] graphically dispalys resistance, compliance and breathing state
- Prompts abnormal compliance and resistance



- Setting range indicator
- Setting summary




- Alarm range indicator
- Active alarm with color coded
- Auto alarm limits


Ease of use

The combination of AMV[™] and IntelliCycle[™] Pro enables the ventilator to make automatic adjustment to ventilator settings, reducing the need for clinicians to make repeated, low level, adjustments so they can direct their focus more effectively on other aspects of patient care.

Thanks to the intuitive UI design of SV300, changing of ventilation modes requires only 2 simple steps. Each function is in logical order so that you wouldn't be lost in complex user manual.







Long-lasting turbine

Powerful built-in turbine provides up to 210 L/min of flow that can fulfill most of the critical patient needs. It does not require any annual maintenance due to its exceptional reliability.



Thoughtful functions built-on latest hardware platform

Future proof

Mindray is committed to protecting your investment by relentlessly developing features that can improve patient care and reduce patient cost. With 25 years experience on development of anesthesia and ventilator products, we proudly present you our new member, SV300 ventilator.

SV300 is built on a highly sophisticate platform, base on state of the art industrial hardware. Your investment is further assure, as we always leave space for future upgrade and comply with the latest stringent medical requirements.

* Accordingly to the latest CE guideline, both the inspiration and expiration valves should be detachable and able to sustain high temperature sterilization in order to minimize the risk of cross contamination.

2 detachable valves

Its design is according to the international standard
*ISO 80601-2-12 & ISO 80601.
We are one of the few manufacturers out there in the market that supports autoclavable inspriration valve. This further reduces the risk of cross contamination.



2 types of EtCO₂ options

SV300 supports both mainstream & sidestream options of EtCO₂ measurement, making it a flexible solution, capable of fitting all of your department needs.



additional functions regarding to patient oxygenation

SV300 provides SpO₂
measurement and O₂ therapy,
both proven to be excellent
tools for patient-weaning
strategies. Also, it allows SpO₂
measurement even when
patient monitor is not around



