## MEDITON MT- 400



## **Electrosurgical Unit**

- ▶ Besides basic surgeries, it also support TURP surgery, under-water cutting, bipolar can work in saline.
- ► Can coonect with endoscope, laparoscope etc
- Monopolar and bipolar, totally 7 working modes
- ► Support vessel sealer (for the vessel < 3mm)
- ► Self-checking system when machine being started
- ► Power self-regualtion
- ► Excellent Electromagnetic Compatibility Function
- Fault alarm and error code display
- ► Alarms for misoperation
- ► Foot switch
- ▶ Power from 250~400 W (produce base on requirement)